

# Micro Reaction Calorimetry:

Newer Applications for the Chemical & Pharmaceutical Industry

Steve Stones • Martyn Ottaway • Thermal Hazard Technology

1 North House • Bond Avenue • Bletchley MK1 1SW • UK

Tel +44 1908 646700 info@thtuk.com



All chemical, physical and biological reactions are accompanied by heat change. These reactions, though sometimes subtle, can be measured using calorimetry. This poster aims to show capabilities of the THT  $\mu$ RC [micro Reaction Calorimeter].

The calorimeter consists of a sample and reference cell designed to accommodate 2ml disposable HPCL glass vials. Optional pressure cells rated to 10bar are also available.

Titration measurements can be made using the automated syringe tower delivering  $\mu$ l injections to the magnetically stirred sample vial.


The in-built peltier allows for experimental temperatures in the range  $-5^{\circ}\text{C}$  to  $150^{\circ}\text{C}$ . Scan rates up to  $2^{\circ}\text{C}/\text{min}$  are also feasible.



## Heat Capacity Measurement

Measurement of heat capacity is an integrated function of the  $\mu$ RC.  $C_p$  is determined by directly measuring the amount of heat required to shift the sample temperature

A small temperature step (usually in the order of  $0.5\text{-}1^{\circ}\text{C}$ ) is applied to the system and the heat is measured. The experiment is then repeated in reverse to verify the measurement. The results from each shift direction are averaged to give the final result. Measurement with empty vials was conducted first to ensure that any differences between the heat capacity of the vials is accounted for.



## Results

The data in the table below show the results from  $C_p$  measurement of pure materials. The values all agree well with the literature data which was obtained from NIST.

	Heat Capacity (Measured)	Heat Capacity (Literature)	Error%
Acetone	2.157	2.130	1.2%
Toluene	1.722	1.704	1.1%
Soya Oil	1.970	1.970	0.0%
NaCl(s)	0.864	0.854	1.2%

The largest error was 1.2% which shows the accuracy of the instrument using sample quantity 1g or less.


This very simple approach makes the THT  $\mu$ RC ideal for rapid (less than 1 hour) measurement of heat capacity of liquids, solids or mixtures.

## Heat Production from Microbes

The sensitivity of the  $\mu$ RC is such that microbial activity can be accurately monitored. Rate of heat production of rumen microbes was successfully studied at  $39^{\circ}\text{C}$  in the calorimeter. Rumen fluid was collected from Jersey cows and the cell suspension (1 mL) added to the sample vial against a reference cell filled with water (1 mL).

## Heat Production from Microbes


The titration syringe was used to deliver 250  $\mu\text{L}$  of 5 mM glucose.




Heat production measured at 1second intervals. Data courtesy of Timothy J. Hackmann, The Ohio State University. USA.

## Heat of reaction under controlled gas flow

The Gas Flow Option (GFO) gives added versatility to the instrument making measurement of the heat of reaction at different flow rates and gas pressures possible.



The option has been used to good effect to study Carbon Capture involving the exothermic reaction that occurs when  $\text{CO}_2$  gas is absorbed by an amine solution. The GFO consists of a flow controller to regulate the flow rate of  $\text{CO}_2$  into the cell containing the amine solution. Weighing the vial before and after the test is used to calculate the rate of  $\text{CO}_2$  uptake.



Methylethanolamine (MEA) was added to water in 30% concentration by weight  $\text{CO}_2$  absorption was monitored at three different flow rates.

## Enthalpy calculated by:

$$\text{Enthalpy of Absorption} = \frac{\text{Heat}}{\text{Number of moles of } \text{CO}_2 \text{ absorbed}}$$

$$\Delta H_{\text{abs}} = \frac{222.4 \text{ J}}{\left(\frac{0.12}{44.01}\right) \text{ mol}} = 81.55 \text{ kJmol}^{-1}$$

Sample	Flow rate of $\text{CO}_2$ feed (ml/min)	$\text{CO}_2$ Absorbed (g)	$\text{CO}_2$ Absorbed (mmol)	Energy released (J)	Enthalpy (kJ/mol)
MEA 0.3 (wt)	0.98	0.12	2.7	222.4	81.55
	1.11	0.12	2.7	226.0	82.87
	0.56	0.11	2.5	209.0	83.60
	0.26	0.12	2.7	224.0	82.13
Average	0.73	0.12	2.7	220.4	82.54


The new Gas Flow option accurately controls  $\text{CO}_2$  dosing and allows direct calculation of  $\text{CO}_2$  loading in the solvent. The ability of the  $\mu$ RC Micro Reaction Calorimeter to quickly carry out heat of absorption and heat capacity measurements using micro litre volumes of reagents makes it ideal for routine carbon capture studies

## Adsorption: $\text{CO}_2$ and $\text{N}_2$ Gas Flow Through Zeolite

The Micro Reaction Calorimeter from THT can perform isothermal enthalpy measurements on the exothermic reactions that occur when gas is absorbed by Zeolites or other adsorbents.

A comparison of the adsorption capacity of an adsorbent with saturated or unsaturated active sites can be made. Using the Zeolite without prior preparation means that it will already have adsorbed atmospheric species and thus its adsorption capacity will be lower than that of regenerated Zeolite. This is demonstrated in the following graph.

In the two tests shown, 100mg of regenerated Zeolite was measured for heat of adsorption in the  $\mu$ RC with a constant  $\text{CO}_2$  flow. The test was then repeated with 100mg of unprepared Zeolite. As expected, the unprepared Zeolite with partially-occupied active sites had a measurably lower adsorption capacity than the regenerated Zeolite.



thermal hazard technology